PROVOX *ActiValve*®

Voice prosthesis for rehabilitation after total laryngectomy

Clinician's
Manual

ATOS

MEDICAL

Leading the way in ENT
















Fig. 16

Emergency Card

Provox ActiValve Lubricant


Fig. 17 a


Fig. 17 b


Fig. 17 c

ENGLISH.....	7
DEUTSCH.....	17
NEDERLANDS.....	27
FRANÇAIS.....	37
ITALIANO.....	47
ESPAÑOL.....	57
PORTUGUES.....	67
SVENSKA.....	77
DANSK.....	87
SUOMI.....	97
HRVATSKI.....	107
繁體中文.....	117
136.....	العربية
Ordering information.....	140

ENGLISH

The Provox ActiValve® Voice Rehabilitation System

Prescriber information

CAUTION: United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

Disclaimer

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Patents and Trademarks

Provox ActiValve is a trademark owned by Atos Medical AB.

Provox ActiValve is protected by US 5.314.470, US 5.976.151, US 7.166.128 and other patents.

This Manual contains information that you may find helpful in considering the use of the Provox ActiValve Voice Rehabilitation System and in assisting you with the use of the device. It is not intended to provide a comprehensive discussion of the risks and benefits of the device, nor is it intended to take the place of a discussion about the risks and benefits of the device between clinician and patient.

The manual, which accompanies this product, may be revised from time to time and must therefore be reviewed prior to each procedure in which the product is used

Stomal airway management

Emergency situations

It is important that emergency personnel knows that you are a neck breather. This will help them determining that oxygen and rescue breathing need to be administered to your stoma and not to your mouth or nose.

We recommend that you and your clinician download the guidelines for rescue breathing from www.atosmedical.com/rescuebreathing

In the package of your Provox ActiValve you will find a bright yellow Emergency card and Provox ActiValve User Cards. We recommend that you keep both the Emergency Card and Provox ActiValve User Card with you at all times

Elective situations

If you need to undergo a procedure that requires intubation (putting a breathing tube in your windpipe), it is very important that the anesthesiologist and the doctor who is performing the procedure are aware that you are a neck breather and that you are using a voice prosthesis. It is important for them to understand that the voice prosthesis should stay in place. If it is removed, fluids from your mouth, food pipe, or stomach may get into your windpipe.

It is very important that the intubation tube is inserted and removed carefully so that the voice prosthesis does not dislodge or come out.

This product and its packaging are latex-free.

Contents

1. Descriptive Information.....	9
1.1 Indications for use.....	9
1.2 Contraindications.....	9
1.3 Description of the device.....	9
1.4 WARNINGS.....	10
1.5 PRECAUTIONS.....	10
1.6 Adverse events.....	10
1.6.1 Prosthesis related.....	10
1.6.2 Puncture related.....	11
2. Instructions for Use.....	13
2.1 Preparation.....	13
2.2 Loading the voice prosthesis into the insertion tool: Figure 1-7 Page 2.....	13
2.3 Anterograde replacement of the prosthesis; Figures 8-15.....	14
2.4 Cleaning and Maintenance.....	14
2.5 Device lifetime.....	15
2.6 Accessories.....	15
2.7 Additional devices.....	15
3. Troubleshooting Information.....	16
3.1 During insertion.....	16
3.2 During use of the ActiValve.....	16
4. Additional Information.....	16
4.1 Warranty.....	16
4.2 Air travel.....	16
4.3 Date of Printing.....	16
4.4 User Assistance Information.....	16

1. Descriptive Information

1.1 Indications for use

The Provox ActiValve is an unsterile indwelling voice prosthesis intended for anterograde insertion in a healed puncture for voice rehabilitation after total laryngectomy. The device is intended for patients who are experiencing early leakage with previous voice prostheses (device life less than 4-8 weeks). The device reduces the need for frequent replacements in a majority of users, but not in all.

1.2 Contraindications

The Provox ActiValve is NOT intended:

- for insertion in a freshly made puncture,
- to be in situ during MRI-examination, (Magnetic Resonance Imaging),
- to be in situ during Radiation Therapy.

1.3 Description of the device

Provox ActiValve voice prostheses are available in one diameter, several lengths (distance between flanges), and three configurations (each having a different given opening force). See graph on inner back cover.

The pressure needed to open the valve flap is primarily controlled by magnets and is different from the pressure needed to open the valve flap in conventional prostheses that do not contain magnets. The magnetic force applies only in the closed and nearly closed position of the valve flap. After the valve flap has opened, the airflow resistance drops. Hence, a higher opening pressure does not mean more air pressure is required during speech. Provox ActiValve Light needs the least pressure to open. Provox ActiValve Strong needs a stronger opening pressure. Provox ActiValve XtraStrong needs the strongest opening pressure. Provox ActiValve XtraStrong is available upon request.

All Provox ActiValve prostheses have the same fluoroplastic material of the valve flap and valve seat. The fluoroplastic material and magnets are radio-opaque and therefore visible on X-ray. The device has no known interference with nuclear or ultrasound imaging techniques (see also Contraindications).

The voice prosthesis is made of medical grade silicone rubber and fluoroplastic, and nickel coated magnets sealed in epoxy glue.

The nickel coated magnets are completely enclosed in the epoxy glue with no detectable release of nickel. The device can therefore also be used in patients with

hypersensitivity to nickel. Care should however be taken not to physically damage the epoxy seal during removal or other manipulation of the device.

Provox ActiValve is manufactured in a cleanroom and supplied NON-sterile, for single use only. It can not be sterilized.

Figure 16, Page 3

A: Inserter

B: Loading Tube

Figure 17, Page 4: Provox ActiValve Voice Prosthesis

- (1) Valve flap, fluoroplastic (radio-opaque); silicone hinge.
- (2) Esophageal flange.
- (3) Tracheal flange.
- (4) Valve seat, fluoroplastic (radio-opaque)
- (5) Interacting magnets in valve seat and in valve flap (radio-opaque).
- (6) Safety strap.
- (7) Direction identification knobs.
- (8) Size information.
- (9) Outer diameter: 7.5 mm (22.5 Fr).

Provox ActiValve package

Each Provox ActiValve comes with the following items:

- 1 single use Insertion Tool (figure 16, A and B) for anterograde insertion,
- 1 bottle (8 ml) Provox ActiValve Lubricant for maintenance,
- 2 reusable Provox Brushes for cleaning of the inside of the prosthesis and for application of Provox ActiValve Lubricant and for inserting the Provox Plug,
- 1 Provox Plug (first aid tool for leakage),
- 2 Provox ActiValve User Cards to be filled in by the clinician,
- 1 Emergency Card,
- 1 Patient's Manual,
- 1 Clinician's Manual.

The brushes, the lubricant, the user cards, the emergency card, the plug and the patient's manual shall be given to the patient.

Provox ActiValve User Card and Emergency Card

Each Provox ActiValve comes with two user cards and one emergency card. The user card provides important information about the Provox ActiValve prosthesis. Please complete the user cards using a water resistant permanent marking pen, provide the cards to the patient.

The emergency card provides important information about emergency airway management. Advise the patient to carry a user card and the emergency card with them at all times in case of emergency, consultation with healthcare providers (e.g. radiology), or during travel. If your patient needs to replace any of the cards, you can order a new one by contacting Atos Medical.

1.4 WARNINGS

When loading the prosthesis in the loading tube always make sure that you fold the esophageal flange **forward** into the loading tube. Incorrect loading can cause immediate valve failure and fluid aspiration.

- Dislodgement or extrusion of the Provox ActiValve voice prosthesis from the TE puncture and subsequent ingestion, aspiration, or tissue damage may occur. For further information about these events, see section Adverse Events.
- Never try to insert the voice prosthesis if its safety strap is damaged or absent. Inserting the prosthesis without the safety strap may increase the risk of aspiration or ingestion of the device.
- Do not modify the Provox ActiValve voice prosthesis or any of its accessories.
- Select the proper prosthesis length. A too tight fit may cause tissue necrosis and extrusion.
- Re-use and re-processing may cause cross-contamination and damage to the device, which could cause patient harm.

1.5 PRECAUTIONS

- Carefully assess patients with bleeding disorders, or patients who are undergoing anticoagulant treatment, for the risk of bleeding or hemorrhage prior to placement or replacement of the prosthesis.
- Aseptic technique should be employed in the management of the voice prosthesis to avoid introducing contaminants into the stoma, TE-puncture or esophagus.
- For internal lubrication of the Provox ActiValve prosthesis (maintenance), only use Provox ActiValve Lubricant. Other lubricants may damage the prosthesis irreparably, leading to material destruction and valve malfunction.

- In patients with excessive biofilm formation (microorganisms and yeasts) the prosthesis, and especially the silicone parts of the prosthesis, may develop heavy biofilm formation. Therefore it is important to instruct the patient to review the Adverse Events section in the manual and to contact their clinician if they experience any of the described events.

1.6 Adverse events

The following are possible risks associated with prosthetic voice rehabilitation:

1.6.1 Prosthesis related

Leakage through the valve – Leakage through the prosthesis may occur because:

- Candida overgrowth near the valve seat and valve flap, or on the silicone hinge may lead to incomplete closure of the valve flap, which causes leakage through the device. This is a normal event in prosthetic voice rehabilitation and an indication to replace the voice prosthesis.
- Excessive Candida overgrowth on the silicone hinge may cause the hinge to break and the flap to be ingested, which causes leakage through the device.
- Stronger negative pressure in the PE-segment occurs during swallowing. To investigate this, trans-stomal inspection of the valve flap of the prosthesis should be performed during swallowing. If the valve opens during swallowing, negative pressure may be the cause of the valve leakage, a Provox ActiValve with a higher magnet strength should be considered.

Candida overgrowth – Candida overgrowth of voice prostheses occurs in almost all patients. In the Provox ActiValve the silicone hinge and prosthesis housing are most prone to be affected by Candida.

Studies suggest that adding certain food items to the laryngectomized patients' diet may reduce the occurrence of yeasts and bacteria in voice prosthesis biofilms. Food items that have been suggested to have such an effect are *buttermilk*^{1,2}, a probiotic milk-drink containing *Lactobacillus Casei Shirota (Yakult)*², *probiotic bacteria Lactobacillus lactis 53* and *Streptococcus thermophilus B*^{3,4} and *Lactobacillus fermentum B54*, *Lactobacillus rhamnosus 744*, *Lactobacillus lactis cremoris SK11*, and *Lactobacillus casei Shirota*⁴, and regular caffeinated soft drinks ("Cola")⁵.

Using the Provox Flush to flush the voice prosthesis with water, blow air through it and imitated coughing may also reduce biofilm formation.⁶

1.6.2 Puncture related

Aspiration of the prosthesis – Accidental aspiration of the voice prosthesis or other components of the voice rehabilitation system may occur. As with any other foreign body, complications from aspiration of a component may cause obstruction or infection. Immediate symptoms may include coughing, wheezing or other abnormal breathing sounds, dyspnea, and respiratory arrest, partial or inadequate air exchange and/or asymmetrical chest movement with respiration. Complications may include pneumonia, atelectasis, bronchitis, lung abscess, bronchopulmonary fistula and asthma. If the patient can breathe, coughing may remove the foreign body. Partial airway obstruction or complete airway obstruction requires immediate intervention for removal of the object.

If aspiration of the device is suspected, a CT scan of the lungs should be performed to confirm aspiration and locate the device. If the CT scan confirms aspiration of the device, the device may be retrieved endoscopically using a non-toothed grasping forceps. The silicone housing of the ActiValve can also be located endoscopically. On a CT scan and during endoscopy, the device may appear as a circular shape with an opening in the middle with an outer diameter of about 14 mm (the flanges of the device), or as a cufflink shape with a shaft length of 4.5, 6, 8, 10 or 12.5 mm, depending on the size of the device. During endoscopy, reflections from the light source on the clear silicone rubber may be seen. Also, in prostheses that have been in situ for some time, white or yellow appearing Candida deposits may be visible on the device.

Ingestion of the prosthesis – Accidental ingestion of the voice prosthesis, or other components of the voice rehabilitation system, may occur. As with any other

foreign body, the symptoms caused by ingestion of the prosthesis or a component of the voice rehabilitation system depends largely on size, location, degree of obstruction (if any) and the length of time it has been present. Ingested components that have remained in the lower esophagus may be removed by esophagoscopy or observed for a short period of time. The object may pass spontaneously into the stomach. Foreign bodies that pass into the stomach usually pass through the intestinal tract. Surgical removal of foreign bodies from the intestinal tract must be considered when bowel obstruction occurs, bleeding is present, perforation occurs or the object fails to pass through the intestinal tract.

Spontaneous passage of the device may be awaited for 4-6 days. The patient should be instructed to observe the stools for the ingested device. If the device does not pass spontaneously, or if there are signs of obstruction (fever, vomiting, abdominal pain) a gastroenterologist should be consulted.

The silicone housing of the ActiValve can be located and retrieved endoscopically. The device may be retrieved by using a non toothed grasping forceps. During endoscopy, the device may appear as a circular shape with an opening in the middle with an outer diameter of about 14 mm (the flanges of the device), or as a cufflink shape with a shaft length of 4.5, 6, 8, 10 or 12.5 mm, depending on the size of the device. Reflections from the light source on the clear silicone rubber may be seen. In prostheses that have been in situ for some time, white or yellow appearing Candida deposits may be visible on the device.

Hemorrhage/Bleeding of the puncture – Slight bleeding from the edges of the TE-puncture may occur during replacement of the prosthesis and generally resolves spontaneously. Patients on anti-coagulant therapy should be carefully evaluated for the risk of hemorrhage prior to placement or replacement of the prosthesis.

Infection and/or edema of the TE-puncture – Infection, granulation formation and/or edema of the puncture may

1. Busscher HJ, Bruinsma G, van Weissenbruch R, Leunisse C, van der Mei HC, Dijk F, Albers FWJ. The effect of buttermilk consumption on biofilm formation on silicone rubber voice prostheses in and artificial throat. *Eur Arch Otolaryngol* 1998; 255: 410-413.

2. Free RH, van der Mei HC, Dijk F, van Weissenbruch R, Busscher HJ, Albers FWJ. Biofilm formation on voice prostheses: influence of dairy products in vitro. *Acta Otolaryngol (Stockh)* 2000; 120: 92-99.

3. Free RH, Busscher HJ, Elving GJ, van der Mei HC, van Weissenbruch R, Albers FWJ. Biofilm formation on voice prosthesis: in vitro influence of probiotics. *Ann Otol Rhinol Laryngol* 2001; 110: 946-951.

4. Van der Mei HC, Free RH, Elving GJ, van Weissenbruch R, Albers FWJ, Busscher HJ. Effect of probiotic bacteria on prevalence of yeasts in oropharyngeal biofilms on silicone rubber voice prostheses in vitro. *J Med Microbiol* 2000; 49: 713-718.

5. Free RH, Elving GJ, van de Mei HC, van Weissenbruch R, Albers FWJ, Busscher HJ. Caffeinated soft drinks reduce bacterial prevalence in voice prosthetic biofilms. *Biofouling* 2000; 120: 69-76.

6. Free RH, van der Mei HC, Elving GJ, van Weissenbruch R, Albers FW, Busscher HJ. Influence of the Provox Flush, blowing, and imitated coughing on voice prosthetic biofilms in vitro. *Acta Otolaryngol (Stockh)* 2003; 123: 547 – 551.

increase the length of the puncture tract. This may cause the prosthesis to be drawn inward and under the tracheal or esophageal mucosa. Inflammation or overgrowth of the esophageal mucosa may also cause the prosthesis to protrude from the puncture. Temporary replacement of the prosthesis by a prosthesis with a longer shaft is then advisable. If standard medical treatment does not resolve the infection, the prosthesis should be removed. In some cases stenting the puncture with a catheter might be considered. If the puncture closes spontaneously repuncture for insertion of a new prosthesis may be required.

Granulation around the puncture – Formation of granulation tissue around the TE-puncture has been reported at an incidence of about 5%. Electrical, chemical, or laser cauterization of the area of granulation may be considered.

Granulation/Hypertrophic scarring around the puncture – Bulging of the tracheal mucosa over the tracheal flange may occur if the prosthesis is relatively short. This excess tissue may be removed by using a laser. Alternatively, a prosthesis with a longer shaft can be used.

Protrusion/extrusion of the prosthesis – Protrusion of the prosthesis and subsequent spontaneous extrusion is sometimes observed during infection of the TE-puncture. Removal of the prosthesis is required to avoid dislodgement into the trachea. The puncture may close spontaneously secondary to the removal of the prosthesis. Repuncture may be necessary for insertion of a new prosthesis.

To reduce the risk of dislodgment/extrusion and its potential sequelae:

- Select the proper prosthesis size.
- Instruct the patient to consult a physician immediately if there are any signs of tissue edema and/or inflammation/infection. Treatment may require antibiotic therapy, temporary insertion of a conventional prosthesis, (e.g. Provox, Provox2 or Provox Vega) with a longer shaft, removal of the prosthesis, surgical closure of the puncture tract, and/or repuncture.
- If laryngectomy tubes or stoma buttons are used, choose a shape that does not exert pressure on the prosthesis and that does not hook the tracheal flange of the prosthesis during insertion and removal of the device. A fenestrated laryngectomy tube (for

speech) should have several small perforations instead of one large opening.

Leakage around the prosthesis – Transient leakage around the prosthesis may occur and may resolve spontaneously. The most common reason is that the prosthesis is too long, which is solved by inserting a shorter prosthesis. If the problem is not solved by inserting the correct length prosthesis, other factors that may affect tissue integrity in the area of the TE puncture (for example gastroesophageal reflux or thyroid function) should be considered, evaluated, and treated. In enlarged punctures with reduced retention strength other conventional treatment methods such as the injection of fillers (e.g. collagen) or temporary removal of the voice prosthesis, should be considered. If the leakage around the voice prosthesis is intractable, more conservative measures, such as surgical closure of the puncture may be necessary.

Tissue damage – If the prosthesis is too short, too long, or is pushed frequently against the esophageal wall by a laryngectomy tube, stoma button, or the patients' finger, damage of the puncture, tracheal and/or esophageal tissues may occur. Inspect the conditions regularly to avoid severe damage.

Instruct your patient to consult a clinician if any of the following occurs:

- Any changes in the appearance of the material of the prosthesis or the way it fits in the puncture,
- leakage during drinking and cleaning the prosthesis does not help,
- speaking becomes difficult (requires higher effort) and cleaning does not help,
- there are any signs of inflammation or tissue change at or near the puncture tract,
- bleeding or overgrowth of tissue around the device,
- persistent pain or discomfort in the region of the Provox ActiValve,
- chronic cough, trouble breathing, or blood in the mucus. These can be signs of a serious health conditions, for example rare lipid pneumonia, or cancer that requires medical attention.

2. Instructions for Use

Caution: Do not use the product if the package is damaged or opened.

If you have little experience with the Provox ActiValve replacement procedure, read the Trouble Shooting section below and practice (optional) on the Artificial Fistula Model Provox (REF 8201). The Artificial Fistula Model consists of a training model with three holes, used to practice insertion of all Provox Voice Protheses.

2.1 Preparation

Choosing the right type of Provox ActiValve

If the present conventional voice protheses valve flap opens during swallowing or heavier breathing, select Provox ActiValve Strong. Otherwise, initially select Provox ActiValve Light.

If an early leakage problem persists with the type of ActiValve chosen, replace the prosthesis with the ActiValve with the next higher opening force, i.e., proceed from Provox ActiValve Light to Strong, or from Provox ActiValve Strong to XtraStrong.

If an upgrade to an ActiValve with a higher opening force is perceived to negatively impact a patient's speech, the patient may return to the use of a device with a lower opening pressure, taking shorter device life into account.

Choosing the right prosthesis type for replacement

If the performance of the current Provox ActiValve was satisfactory, the same type should be used for replacement. If the device lifetime of the current prosthesis was too short, or if speaking was perceived as too difficult choose a new prosthesis as described above.

Choosing the right length of the prosthesis

To select the correct length, you may use the current prosthesis as its own measuring device. If there is too much (i.e. 3 mm / ~0.12 inches, or more) space between the Tracheal Flange of the current prosthesis and the mucosal wall, a shorter prosthesis should be used for replacement.

If the prosthesis sits too tight, a longer prosthesis should be used. Deformation of the tracheal flange may indicate that the prosthesis is too short. Use of a longer prosthesis should be considered if the tracheal mucosa shows a tendency to become hypertrophic and to bulge over the tracheal flange.

2.2 Loading the voice prosthesis into the insertion tool: Figure 1-7 Page 2

The following method is described for right-handed use; for left-handed use, references to handedness should be reversed.

See Figure (1) The insertion tool consists of an inserter (A) and a loading tube (B). The inserter has a key hole opening for securing the prosthesis. The loading tube has a tapered end for anterograde insertion into the TE-puncture. The loading tube has been pre-lubricated to facilitate passage of the voice prosthesis through the tube. Insert the safety strap of the prosthesis into the key hole opening (2) of the inserter, fixate it by pulling firmly in the direction of the arrow (3) and place the prosthesis with its tracheal opening on top of the inserter (4).

(5) Squeeze the esophageal flange of the voice prosthesis forward, between the thumb and index finger of your right hand, in the direction that will be toward the esophagus during device placement. Hold the loading tube with your left hand. Insert ("hook") the squeezed flange at an angle of about 45° into the slit of the loading tube.

WARNING: If the esophageal flange is not folded in the direction that will be toward the esophagus during device placement it could damage the device e.g., by dislodging the blue valve seat or tearing off the valve flap. If the valve is damaged do not use the device. There is a risk of aspiration.

(6) Push the inserter forward for about 1 cm (~0.4 inches). This holds the voice prosthesis with most of the esophageal flange folded forward (flanges bent into the loading tube).

(7) Press the protruding part of the flange down with your left thumb in to the loading tube. The esophageal flange should now be completely folded forward into the loading tube, pointing towards the front tip of the loading tube. Push the inserter forward until mark line 1 is level with the rim of the loading tube.

The Provox ActiValve prosthesis is now prepared for placement.

If the esophageal flange is not completely folded forward into the loading tube and it is hard to push the prosthesis forward within the tube, interrupt the insertion and pull the prosthesis backwards out of the loading tube. Reload and repeat the procedure. Additional water-soluble lubrication may be needed at the third attempt.

2.3 Anterograde replacement of the prosthesis; Figures 8-15

Remove the old prosthesis with the help of non-toothed hemostats

See Figure (8): Insert the tapered side of the loading tube into the TE-puncture, and probe the puncture pushing forward, until the tip of the loading tube is completely inside the lumen of the esophagus and the posterior esophageal wall is reached.

(9) Hold the loading tube with your left hand in this position. Hold the inserter with your right hand, using your thumb as a natural stop at mark line 2. Push the inserter (right hand) forward until mark line 2 is level with the rim of the loading tube.

CAUTION: Do not “inject” the prosthesis using the inserter like a syringe. Do not hold the inserter like a pen. Both methods will lead to incorrect placement.

(10) The forward-folded esophageal flange should now unfold in the lumen of the esophagus. Check this by gently pulling the insertion tool (loading tube and inserter together) backwards until the anterior esophageal wall is reached. Hold the inserter in this position.

(11) Pull the loading tube backwards while keeping the inserter in the same position and remove the tube over the inserter. As a result, the tracheal flange should unfold. Check to see if the prosthesis is in the proper position by rotating and exerting slight traction on the prosthesis. If the tracheal flange does not unfold automatically, pull lightly and turn the inserter with the prosthesis around its axis several times, or use a non-toothed hemostat to pull gently on the tracheal flange and rotate this in position. The knobs on the prolonged tracheal flange indicate in which direction the safety strap should be pulled. The smooth side of the tracheal flange should face the tracheal mucosa and the knobs should face toward the lumen.

(12) When the prosthesis is correctly in place, cut off the safety strap on the outer edge of the tracheal flange. The end of the tracheal flange should point downwards into the trachea.

(13, 14) Lubricate the inner lumen of the prosthesis with Provox ActiValve Lubricant, using one drop on a Provox Brush. While doing so, it is advisable to explain and demonstrate maintenance procedures to the patient.

(15) The Provox ActiValve voice prosthesis is ready for use.

14

If the insertion is not successful, it can be repeated with the same tools and devices. If, however, the procedure has to be repeated more than twice, additional water-soluble lubrication of the loading tube shall be used.

2.4 Cleaning and Maintenance

Caution: Only use genuine Provox accessories that are intended for use with Provox ActiValve for the handling and cleaning of your prosthesis.

Instruct your patient on how to clean the device and to avoid all other kinds of manipulation of the prosthesis than the cleaning and maintenance described below.

The patient should clean the prosthesis at least twice a day and after each meal with the Provox Brush by inserting the brush into the prosthesis and moving it gently back and forth with a twisting movement. After the brush is removed it should be wiped with a gauze. The procedure can be repeated as often as necessary. For details and how to clean the brush see Instructions for Use that come with the Provox Brush.

In addition to using the Provox Brush, the patient may also clean the Provox ActiValve voice prosthesis with the Provox Flush. The Provox Flush can be used with drinking water or air. For details and how to clean the Provox Flush see Instructions for Use that come with the Provox Flush.

Cleaning and disinfection of the accessories

The accessories should be cleaned after each use and disinfected at least once a day according to the Instructions for Use.

Inform the patient if the accessories are visibly soiled or have air dried when there is someone in the home with a respiratory infection, or have been otherwise contaminated by contact with a pet, dropping on the floor or other gross contamination, the accessories should be cleaned and disinfected just prior to use.

During hospitalization there is an increased risk of device contamination and patient infection. Therefore it is important to clean and disinfect the accessories immediately after use and again just before use during hospitalization, and rinse using sterile water, rather than drinking water.

Lubrication

Once a day, after cleaning in the morning, the inner lumen of the Provox ActiValve prosthesis should be

lubricated with one drop of Provox ActiValve Lubricant using a Provox Brush. Lubricating the inner lumen of the prosthesis helps prevent occasional stickiness of the valve that might otherwise occur, typically after sleep or a prolonged time of not speaking.

The patient should apply only a single drop of the Provox ActiValve Lubricant by using a different Provox Brush. The lubricant bottle releases the lubricant very slowly, one drop at a time. It is important to be patient when squeezing the Lubricant bottle. The liquid is thick and will come out slowly. Instruct the patient to not put more than one drop on the brush and avoiding touching the bristles of the brush with the Provox ActiValve Lubricant bottle when dispensing the drop.

After use, wipe the brush thoroughly with a gauze. For how to clean the brush see Instructions for Use that come with the Provox Brush.

The Provox ActiValve Lubricant can be used until the expiration date or 8 months after opening the bottle, whichever occurs first. Inform the patient to register the date of opening on the bottle using a water resistant permanent marking pen.

Note: Incidental ingestion of small amounts of Provox ActiValve Lubricant is harmless. Incidental aspiration of small amounts of Provox ActiValve Lubricant is also harmless, but may cause coughing.

No health problems have been reported regarding repeated aspiration of small amounts of silicone oil, but it is still better to be cautious and avoid it. Use one drop of oil each time and make sure that any superfluous oil is removed before using the brush (no dripping brush). Avoid touching the stomal mucosa with the brush during application of the Provox ActiValve Lubricant. Remove any excess oil from the brush after use by wiping it with a gauze.

2.5 Device lifetime

The prosthesis is not a permanent implant, and requires periodic replacement. Depending on individual biological circumstances the device life varies, and it is not possible to predict the integrity of the device over a longer period of time. The prosthesis and especially the silicone material of the device will be affected by e.g. bacteria and Candida, and the structural integrity of the device will eventually deteriorate.

Indications for replacement of the Provox ActiValve voice prosthesis include leakage through the valve, blockage of the prosthesis, bacterial and candidal overgrowth leading to degradation of the materials and/or excessive pressure needed to obtain speech.

Other reasons for earlier replacement could include medical indications such as problems with the puncture tract (see Adverse Events section).

Compatibility with antifungal medicine

In most cases treatment with antifungal agents should not be indicated but may be considered as a preventive measure if excessive candida overgrowth occurs on the prosthesis.

An unknown variety of chemical substances may influence the material properties of the device. Therefore, the introduction of antifungal medicines or drugs directly to or within close proximity of the voice prosthesis should be carefully assessed.

Laboratory testing shows no negative influence on the function of the Provox ActiValve voice prosthesis and components when using the following antifungal medications: Nystatin, Fluconazole and Miconazol.

2.6 Accessories

CAUTION: Only use genuine Provox accessories that are intended for use with Provox ActiValve. Other accessories may cause patient harm or product malfunction.

- Accessories to facilitate insertion of the prosthesis (for clinician's use)
Provox Dilator: Helps to convert from a prosthesis with a smaller shaft diameter to a Provox voice prostheses of larger diameter.
- Accessories to maintain the device function (for patient's use)
Provox Brush/Provox Brush XL is used by the patient to clean the inside of the prosthesis.
Provox Flush: An additional cleaning device for patient use that allows flushing of the prosthesis.
Provox Plug: A tool for patient's use that temporarily blocks leakage through the prosthesis.
Provox ActiValve Lubricant: Serves as lubricant for the valve surface inside the Provox ActiValve prosthesis.

2.7 Additional devices

- **Provox HME:** Heat- and Moisture Exchanger. Partially restores lost nasal functions (heating, humidifying, breathing resistance).
- **Provox FreeHands HME:** A Heat and Moisture Exchanger combined with an automatic tracheostoma valve, designed to allow hands-free speech.

- **Provox Adhesives:** Different types of adhesives to attach HME's.
- **Provox LaryTube:** A silicone laryngectomy tube intended to support shrinking stomas while retaining other rehabilitation devices of the Provox System.
- **Provox LaryButton:** A soft, self-retaining silicone stoma button intended to support shrinking stomas while retaining other rehabilitation devices of the Provox System.
- **Artificial Fistula Model Provox:** A training model that consists of a plate with a holder which has three holes, used to practice insertion of all Provox Voice Prostheses.

3. Troubleshooting Information

3.1 During insertion

Inserting the prosthesis completely into the esophagus (overshooting).

In case of overshooting (the whole prosthesis is developed in the esophagus), leave the safety strap attached to the inserter and use a non-toothed hemostat to pull the prosthesis into position. The identification knobs on the prolonged tracheal flange must always face the open lumen of the trachea, never the mucosal membrane. Slight rotation might be helpful to fully position the prosthesis. Cut the safety strap only after the correct fit is obtained.

Forgotten to stop at mark "2"

If this happens, see the section above on "overshooting."

3.2 During use of the ActiValve

- **Symptom:** Speaking is sometimes not possible (valve does not open).
Probable cause: The valve lid sticks to the valve seat, preventing the valve from opening. This becomes more likely the longer the valve has not been in use (e.g. overnight).
What to do: Clean the prosthesis and lubricate it (see section Cleaning and Maintenance). If cleaning and lubricating does not help, the prosthesis needs to be exchanged.
- **Symptom:** Coughing triggered by drinking or eating while the prosthesis is in place.
Probable cause: Leakage through or around the voice prosthesis.
What to do: Clean the prosthesis. There may be food remnants preventing the valve from proper closure. **If cleaning does not help:** Close the prosthesis with a Provox Plug. If the leakage stops, the prosthesis needs to be replaced.

If the problem persists, even with the Provox

Plug in place: There is leakage around the prosthesis. See section Adverse events - Leakage around the prosthesis.

- **Symptom:** Speaking becomes increasingly difficult.
Probable cause: Valve obstruction, tissue swelling in the throat.
What to do: Clean the prosthesis. There may be food remnants increasing airflow resistance. If cleaning does not help the prosthesis needs to be exchanged.

4. Additional Information

4.1 Warranty

See disclaimer section.

4.2 Air travel

This device contains 2 magnets; Size 1-1.5 mm. weight 0.15 – 0.2 gram. The magnets are unlikely to set off metal detectors such as those used in airports but in order to avoid incidents it is advisable to inform the security personnel in advance of the voice prosthesis and show them the user card.

4.3 Date of Printing

See version number on the back cover of this manual.

4.4 User Assistance Information

For additional help or information please see back cover of this manual for contact information.

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Atos Medical Inc.