



## 4801N Epiflex

<b>DE</b> Gebrauchsanweisung .....	3
<b>EN</b> Instructions for Use .....	8
<b>FR</b> Instructions d'utilisation .....	12
<b>IT</b> Istruzioni per l'uso .....	17
<b>ES</b> Instrucciones de uso .....	22
<b>NL</b> Gebruiksaanwijzing .....	27

	Größe / Size	Ellenbogenumfang / Elbow circumference	
		cm	inch
	S	24–25	9.51–9.8
	M	26–27	10.2–10.6
	L	28–30	11.0–11.8
	XL	31–33	12.2–13.0

<b>Material</b>	Kunststoffspange / plastic support: PC, PE, SBR (Neopren / neoprene)
	Gewebe / Fabric: PE, EL

# 1 Vorwort

## INFORMATION

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### INFORMATION

Date of the last update: 2012-11-14

- ▶ Please read this document carefully.
- ▶ Follow the safety instructions.

These Instructions for Use provide you with important information on the handling and use of the 4801N Epiflex epicondylitis orthosis.

## 2 Application

### 2.1 Intended Use

The orthosis is to be used **solely** as an epicondylitis orthosis and is intended **exclusively** for contact with healthy skin.

The orthosis must be used according to the indications and following the instructions of a physician.

### 2.2 Indications

- Epicondylitis humeri radialis (tennis elbow)
- Epicondylitis humeri ulnaris (golfer's elbow)

Indications must be determined by the physician.

### 2.3 Contraindications

#### 2.3.1 Absolute Contraindications

Not known.

#### 2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases / injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the support / orthosis will be applied; sensory and circulatory disorders in the body region being treated and in the area of the forearm and hand.

## 2.4 Effects

The orthosis with a semi-rigid, moulded plastic support has an analgesic effect because it relieves strain on the muscle origins of the forearm muscles (depending on the positioning of the orthosis, the muscle origins of the hand and finger extensors or of the hand and finger flexors). It allows daily activities to be completed, notwithstanding elbow complaints.

The orthosis is universal for the left or right side.

## 3 Safety

### 3.1 Explanation of Warning Symbols

Warning regarding possible risks of accident or injury

Warning regarding possible technical damage

### 3.2 General Safety Instructions

#### **Using the product without following the instructions for use**

Deteriorating health condition and damage to the product because of failure to follow the safety information

- ▶ Observe the safety information in these instructions for use.
- ▶ Make the patient aware of all safety instructions marked with "**Inform the patient**".

#### **Reuse on other patients and improper cleaning**

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Only use the product for a single patient.
- ▶ Clean the product regularly.
- ▶ **Inform the patient.**

**⚠ CAUTION**

**Use of a product with neoprene**

Allergic reactions due to neoprene or heat allergies

- ▶ Do not wear products with neoprene for more than 3 – 4 hours without interruption (unless explicitly instructed otherwise by the doctor).
- ▶ Do **not** use the product in case of a known neoprene or heat allergy.
- ▶ **Inform the patient.**

**⚠ CAUTION**

**Contact with heat, embers or fire**

Injuries (e.g. burns) caused by melting of the material

- ▶ Keep the product away from open flame, embers and other sources of heat.
- ▶ **Inform the patient.**

**NOTICE**

**Contact with salves, lotions or other products that contain oils or acids**

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to salves, lotions or other products that contain oils or acids.
- ▶ **Inform the patient.**

## 4 Handling

**INFORMATION**

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.

- Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

## 4.1 Selecting the Size

- 1) Determine the elbow circumference.
- 2) Determine the orthosis size (see size table).

## 4.2 Fitting and Application

### CAUTION

#### **Incorrect or excessively tight application**

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- Ensure that the orthosis is applied properly and fits correctly.
- **Inform the patient.**

- 1) Open the hook-and-loop closure. Do not fully pull the strap out of the loop in this process.
- 2) Pull the plastic support over the forearm. Ensure that the large plate is on the outside (see Fig. 1, Pos. 1).
- 3) Position the plastic support **approx. 5 cm** below the elbow joint (see Fig. 2).
- 4) Tighten the closure strap so that a finger still fits between the arm and closure strap (see Fig. 1, Pos. 2).
- 5) Close the hook-and-loop closure.

## 4.3 Cleaning

Clean the orthosis regularly.

- 1) Close all hook and loop closures.
- 2) **Recommendation:** Use a laundry bag or net.
- 3) Wash the product in 30 °C warm water with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat (e.g. from the sun, an oven or radiator).

## 5 Disposal

Dispose of the product in accordance with national regulations.

## 6 Legal Information

### 6.1 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

### 6.2 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH