8360 Genu Carezza Patella Stabilizer

DE	Gebrauchsanweisung	3
EN	Instructions for use	7
FR	Instructions d'utilisation	11
IT	Istruzioni per l'uso	15
ES	Instrucciones de uso	19
NL	Gebruiksaanwijzing	24
SV	Bruksanvisning	28
DA	Brugsanvisning	31
NO	Bruksanvisning	35
FI	Käyttöohje	39
PL	Instrukcja użytkowania	43
CS	Návod k použití	47
HR	Upute za uporabu	51
BG	Инструкция за употреба	55
TR	Kullanma talimatı	59
RU	Руководство по применению	63
ZH	使用说明书	67

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INFORMATION

Last update: 2015-02-27

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

These instructions for use provide important information on the fitting and application of the 8360 Genu Carezza patella stabilizer knee support.

2 Application

2.1 Indications for use

The support is intended **exclusively** for fittings of the lower limbs and **exclusively** for contact with intact skin.

The support must be used according to the indication.

2.2 Indications

- Patellofemoral pain syndrome (e.g. chondromalacia patellae, chondropathia patellae, retropatellar osteoarthritis)
- Patellofemoral maltracking
- Following patellar dislocation/subluxation
- After lateral release operations

The indication must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

None known.

2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening

and hyperthermia of the fitted limb/body area; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the medical device will be applied; sensory and circulatory disorders in the leg area.

2.4 Mechanism of Action

The support stabilizes the kneecap laterally and relieves pain.

3 Safety

3.1 Explanation of Warning Symbols

△ CAUTION	Warning regarding possible risks of accident or injury.
NOTICE	Warning regarding possible technical damage.

3.2 General Safety Instructions

△ CAUTION
Contact with heat, embers, or fire
Injuries (e.g. burns) caused by melting of the material
► Keep the product away from open flame, embers, and other sources of heat.

△ CAUTION
Reuse on other patients and insufficient cleaning
Skin irritation, formation of eczema, or infections due to contamination with germs
► Use the product only for a single patient. ► Clean the product regularly.

NOTICE
Improper use and changes
Change in or loss of functionality as well as damage to the product
► Use this product with care and only for its intended purpose. ► Do not make any improper changes to the product.

NOTICE

Contact with ointments, lotions, or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- Do not expose the product to ointments, lotions, or other products that contain oils or acids.

4 Handling

INFORMATION

- The daily duration of use and period of application are generally determined by the physician.
- The initial fitting and application of the product must be carried out by qualified personnel.
- Instruct the patient in the handling and care of the product.
- Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Selecting the Size

- 1) Measure the circumference of the thigh 15 cm above the center of the patella.
- 2) Measure the circumference of the lower leg 15 cm below the center of the patella.
- 3) Determine the size of the support (see size chart).

4.2 Fitting and Application

CAUTION

Incorrect or excessively tight application

Skin irritation due to overheating, local pressure points due to tight application to the body

- Ensure that the support is applied properly and fits correctly.
- Do not use the product in case of a known heat allergy.
- When in doubt, do not continue to use the product if skin irritation occurs.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Instruct the patient to inspect the product for proper functionality, wear, and damage prior to each use.
- ▶ Advise the patient that the product may not be used if the product or part of the product exhibits signs of wear (e.g. cracks, deformation, poor fit) or damage.

- > The patient is sitting on the front edge of a chair.
- > Open all hook-and-loop closures on the support.
- 1) Hold the support by the top edge on both sides and pull it over the knee until the patella cutout is centered on the patella.
- 2) Fasten the straps first below and then above the patella (see fig. 1).

4.3 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the support regularly:

- 1) Fasten all hook-and-loop closures.
- 2) Recommendation: Use a laundry bag or net.
- 3) Wash the support in warm water at **40 °C** (104 °F) with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

5 Disposal

Dispose of the product in accordance with national regulations.

6 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

6.1 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorized modification of the product.

6.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The manufacturer therefore drew up the declaration of conformity on its own responsibility in accordance with Annex VII of the directive.

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Otto Bock HealthCare GmbH