
PROVOX[®] Vega[™] XtraSeal[™]

Voice Prosthesis with SmartInserter[™]
Clinician's Manual

**Rx
ONLY**

CE 0413

STERILE EO


Atos

Atos Medical Your voice

Illustrations


1. Provox® Vega™ XtraSeal™ voice prosthesis

- a) Size information (shaft diameter and length between flanges)
- b) Prosthesis Hood
- c) Esophageal Flange
- d) Tracheal Flange
- e) Safety Strap
- f) Valve Flap
- g) Prosthesis Shaft
- h) Radio-opaque fluoroplastic Valve Seat
- i) Additional enlarged Esophageal Flange


2. SmartInserter™ System

- a) Guide
- b) Connection Side
- c) Folding Device
- d) Distal tip
- e) Inserter Pin
- f) Voice prosthesis
- g) Loading Tube
- h) Step marker
- i) Tip


3. Over-shooting


3.1


3.2

3.3


3.4


3. Over-shooting


3.5


3.6

3.7


3.8


3. Over-shooting

3.9


3.10


4. Re-loading


4.1


4.2


4.3


4.4

4.5

4.6

Prescription information

CAUTION: United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

Disclaimer

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Patents and trademarks

Provox® is a registered trademark owned by Atos Medical AB, Sweden. Provox® Vega™ and Provox® Vega™ XtraSeal™ are protected by US patent no 8.332.999 as well as other patents, patents pending and design patents.

Unregistered copy

Contents

ENGLISH	10
DEUTSCH	16
NEDERLANDS.....	22
FRANÇAIS	28
ITALIANO.....	34
ESPAÑOL	40
PORTUGUÊS.....	46
SVENSKA	52
DANSK.....	57
NORSK.....	63
SUOMI	69
ČESKY	75
POLSKI.....	81
HRVATSKI	88
TÜRKÇE	94
РУССКИЙ	100
日本語.....	107
한국어	113
繁體中文.....	118
עברית.....	129
العربية.....	136
ORDERING INFORMATION.....	137

Provox® Vega™ XtraSeal™

The Instructions for Use, which accompanies this product, may be revised from time to time and must therefore be reviewed prior to each procedure in which the product is used.

Contents

1. Descriptive information.....	10
1.1 Intended use	10
1.2 Description of the device	10
1.3 CONTRAINDICATIONS	10
1.4 WARNINGS	10
1.5 PRECAUTIONS.....	11
2. Instructions for use.....	11
2.1 Preparation	11
2.2 Anterograde replacement procedure using overshooting	11
2.3 Reloading the Provox SmartInserter	12
3. Important patient information	12
3.1 General information	12
3.2 Prosthesis maintenance.....	12
3.3 Accessories.....	13
4. Adverse events/Trouble shooting information.....	13
5. Additional information	14
5.1 Compatibility with MRI, X-ray and radiation therapy.....	14
5.2 Device lifetime.....	15
5.3 User assistance information	15

1. Descriptive information

1.1 Intended use

The Provox Vega Voice Prosthesis is a sterile single use indwelling voice prosthesis intended for voice rehabilitation after surgical removal of the larynx (laryngectomy). Cleaning of the voice prosthesis is performed by the Patient while it remains in situ.

The Provox SmartInserter is a sterile single use device intended for anterograde replacement of the Provox Vega Voice Prosthesis. This replacement procedure is carried out by a medical doctor or a trained medical professional in accordance with local or national guidelines.

The Provox SmartInserter is not intended to be used for insertion of a voice prosthesis in a freshly made puncture.

1.2 Description of the device

General

The Provox Vega is a one-way valve (prosthesis) that keeps a TE-puncture open for speech, while reducing the risk of fluids and food entering the trachea. The Provox Vega voice prosthesis is not a permanent implant, and needs periodic replacement. The prosthesis (Fig. 1) is available in different diameters and several lengths. Provox Vega XtraSeal has an additional enlarged esophageal flange that is intended to solve problems with leakage around the voice prosthesis. The device is made of medical grade silicone rubber and fluoroplastic.

Provox Vega XtraSeal package

The Provox Vega XtraSeal package contains the following items:

- 1 Provox Vega XtraSeal voice prosthesis, pre-loaded in a single-use SmartInserter, sterile (Fig. 2)
- 1 Provox Brush of a size corresponding to the voice prosthesis, non-sterile
- 1 Provox Vega XtraSeal Clinician's manual
- 1 Provox Vega Patient's manual
- 1 Provox Brush Instructions for Use

1.3 CONTRAINDICATIONS

There are no known contraindications for use or replacement of a Provox Vega voice prosthesis among patients already using prosthetic voice rehabilitation.

1.4 WARNINGS

Dislodgement or extrusion of a Provox Vega voice prosthesis from the TE puncture and subsequent ingestion, aspiration or tissue damage may occur. For further information about these events, see section **Adverse Events/Trouble Shooting Information** below.

To reduce the risk of dislodgment/extrusion and its potential sequelae:

- Select the proper prosthesis size. A tight fit may cause tissue necrosis and extrusion.
- Instruct the patient to consult a physician immediately if there are any signs of tissue edema and/or signs of inflammation/infection.
- If used, choose laryngectomy tubes or stoma buttons with a suitable shape that do not exert pressure on the prosthesis or catch onto the tracheal flange of the prosthesis during insertion and removal of the laryngectomy tube or stoma button.
- Instruct the patient to use only genuine Provox accessories of corresponding size (Brush, Flush, Plug) for maintenance, and to avoid all other kinds of manipulation.
- Re-use and re-processing may cause cross-contamination and damage to the device, which could cause patient harm.

1.5 PRECAUTIONS

- Carefully assess any patient with bleeding disorders or who is undergoing anticoagulant treatment for the risk of bleeding or hemorrhage prior to placement or replacement of the prosthesis.
- Inspect the package before use. If damaged or opened, do not use the product.
- Always use aseptic technique when handling the prosthesis to reduce infection risk.
- Make sure that any external or internal stoma attachment devices used do not exert pressure or hook the flange of the prosthesis (e.g., HME base plates, laryngectomy tubes, or stoma buttons). This may lead to severe tissue damage and/or accidental ingestion of the prosthesis.

2. Instructions for use

2.1 Preparation

Choosing the right size of the replacement prosthesis

Be sure to use a Provox Vega XtraSeal voice prosthesis of the proper shaft diameter and length. Provox Vega XtraSeal is available in different length/diameter combinations.

- **Selecting shaft diameter**

The Clinician should determine the proper diameter of the prosthesis appropriate for the patient.

- If the selected diameter is larger than the previous prosthesis the tract must be dilated, using the Provox Dilator, appropriate to the diameter of the prosthesis being inserted.
- If a prosthesis with a smaller shaft diameter is inserted, observe and ensure that the puncture shrinks to the appropriate diameter.

- **Selecting shaft length**

To select the correct length, you may use the current prosthesis as its own measuring device. If there is too much (i.e. 3 mm / ~0.12 inches, or more) space between the Tracheal Flange of the old prosthesis and the mucosal wall, a shorter prosthesis should be used. If the prosthesis sits too tight, a longer prosthesis should be used.

Note: The shaft of Provox Vega XtraSeal is ca 1 mm shorter than the size indicated due to the enlarged esophageal flange.

2.2 Anterograde replacement procedure using overshooting

(Fig. 3.1 – 3.10)

CAUTION: Provox Vega XtraSeal with enlarged esophageal flange must be placed using the overshooting technique to ensure the additional enlarged esophageal flange and the esophageal flange both deploy inside the esophageal lumen. Overshooting means that the whole prosthesis is deployed in the esophagus and then retracted to the intended position.

Activating System for replacement using overshooting

1. With one hand holding the loading tube, use the other hand to push the inserter pin partially into the loading tube, until the blue ring of the prosthesis has passed the step marker of the loading tube (Fig. 3.1).
2. Squeeze the Guide between thumb and index finger (Fig. 3.2). Push it further into the loading tube (Fig. 3.3). The Provox Vega XtraSeal voice prosthesis is now ready for overshoot-insertion.

Inserting the new prosthesis into the puncture - Overshoot

1. **Probing the TE-puncture**

Hold the Provox SmartInserter by the Loading Tube. Gently insert the Tip of the Loading Tube. Proceed with care if you encounter resistance during insertion of the Loading Tube into the puncture. In this case, dilatation and/or lubrication of the puncture can ease the insertion

2. **“Injecting” the prosthesis**

Hold the loading tube stable in this position and use the other hand to advance the Inserter Pin until completely inserted into the SmartInserter (Fig. 3.4). As the Guide is being pushed further into the Loading Tube, the Inserter Pin will advance further and release the Provox Vega XtraSeal voice prosthesis completely within the esophagus.

3. **Pulling the prosthesis back**

Pull the Loading Tube out from the puncture (Fig. 3.5). Pull the Provox Vega XtraSeal voice prosthesis back into the puncture by grasping and pulling the tracheal flange of the prosthesis with a non-toothed hemostat (Fig. 3.6).

4. **Checking the proper fit**

Check to see if the prosthesis and the flanges are in the proper position by rotating and moving the prosthesis back and forth in the puncture. Apply slight traction on the tracheal flange to ensure proper retention of the voice prosthesis (Fig 3.7).

Caution: Confirm that both esophageal flanges have deployed entirely in the lumen of the esophagus. No part of the enlarged esophageal flange should be visible alongside the shaft of the voice prosthesis and when rotated the voice prosthesis should move freely. Rotating the voice prosthesis while at the same time slightly pushing it towards the esophagus can help to unfold the enlarged flange completely. In case of uncertainty, use flexible endoscopy to confirm proper placement inside the esophagus.

5. Test proper function

Test proper function of the prosthesis by asking the patient to speak (Fig. 3.8) and by observing the prosthesis while the patient drinks water (Fig. 3.9).

6. Cut the Safety Strap with scissors

When the prosthesis is confirmed to be in place, position the oval end of the tracheal flange pointing downwards and use a pair of scissors to cut off the Safety Strap on the outer edge of the flange (Fig. 3.10). The voice prosthesis is now ready for use.

2.3 Reloading the Provox SmartInserter

(Fig. 4.1 – 4.6)

If the first attempt to insert the Provox Vega XtraSeal voice prosthesis into the TE-puncture was unsuccessful, it is possible to reload the Provox Vega XtraSeal into the SmartInserter.

CAUTION: Do not reload the SmartInserter more than 2 times.

If the Provox Vega XtraSeal voice prosthesis shows any signs of damage, do not use the prosthesis.

1. Make sure the prosthesis is still safely attached to the Inserter Pin. If not, secure the Safety Strap on the Inserter Pin by leading the Safety Strap through the Attachment Hole in the Inserter Pin (Fig. 4.1).
2. Put the Provox Vega XtraSeal voice prosthesis with the tracheal side down on top of the Inserter Pin (Fig. 4.2). Use your thumb to ensure that the voice prosthesis is properly positioned.
3. Check that the Folding Device and Guide are still assembled in the upper Guide position. If not, re-assemble the Folding Device into Guide (Fig. 4.3). Do not push the Folding Device in the bottom position of the Guide.
4. Gently pull the Inserter Pin with the attached prosthesis into the Folding Device until you can clearly hear it click in place (Fig. 4.4). The Esophageal Flange of the prosthesis is now in level with the Folding Device.
5. The enlarged E-flange must be manually positioned below the curved tips of the folding device (Fig. 4.5).
NOTE: If this step is omitted, preloading the voice prosthesis will require excessive force and may damage the voice prosthesis.
6. Attach the Loading Tube to the top of the Folding Device / Guide (Fig. 4.6) and let the two parts snap together.

The Provox Vega XtraSeal voice prosthesis is now ready for activation and can be used as described in section “**Anterograde Replacement Procedure using overshooting**”.

3. Important patient information

3.1 General information

Ensure that the patient understands to contact their clinician if:

- There is leakage through or around the prosthesis (coughing and/or change of mucus color).
- Speaking becomes difficult (higher effort and/or voice sounds more strained).
- There are any signs of inflammation or tissue changes at the puncture site or the stoma region (pain, redness, heat, swelling, traces of blood on the brush after brushing).

Also inform the patient that:

- After a prosthesis replacement traces of blood may be found in the sputum. This may come from granulation tissue on the edges of the TE puncture.
- Occasionally, mild leakage through or around the prosthesis may occur in the first weeks after insertion of a new prosthesis. This often resolves spontaneously and does not require immediate replacement of the prosthesis.
- If leakage spontaneously occurs through the prosthesis a temporary blockage can be done with a Provox Vega Plug.
- Speech Training Sessions with a trained clinician are advised in order to acquire optimal voice sound, fluent speech, and optimal intelligibility.

3.2 Prosthesis maintenance

Instruct the patient when and how to clean the Provox Vega XtraSeal voice prosthesis. Cleaning the prosthesis at least twice a day can help prolong the device life.

CAUTION: Only use genuine Provox accessories that are intended for use with Provox Vega when cleaning the prosthesis.

- Brushing the prosthesis with the Provox Brush twice a day will help to remove mucus and food remnants from the prosthesis.
- Flushing the prosthesis with the Provox Flush also helps to clean debris and mucus from the prosthesis, which helps increase the life of the device.

Note: The Provox Flush is intended for use only by patients who, as assessed by the clinician who prescribes the device, have demonstrated ability to understand and consistently follow the instructions for use without clinician supervision.

- Some dietary measures, like the daily intake of yogurt or butter milk containing lactobacilli, are considered to be helpful against excessive Candida growth.

For detailed information on how to clean a Provox Vega, please see Instructions for each accessory.

3.3 Accessories

CAUTION: Only use genuine Provox accessories that are intended for use with Provox Vega. Other accessories may cause patient harm or product malfunction.

Accessories to facilitate insertion of the prosthesis (for clinician's use)

- **Provox Dilator:** Helps to convert from a prosthesis with a smaller shaft diameter to a Provox voice prostheses of larger diameter.
- **Provox GuideWire:** for placement of a Provox, Provox2 or a Provox Vega voice prosthesis in a freshly made TE puncture, or for retrograde replacement of a Provox indwelling voice prosthesis.

Accessories to maintain the device function (for patient's use)

- **Provox Brush / Provox Brush XL** is used by the patient to clean the inside of the prosthesis.
- **Provox Flush** is an additional cleaning device that allows flushing of the prosthesis.
- **Provox Vega Plug** is a First Aid tool for patient's use that temporarily blocks leakage through the prosthesis.

Other supplementary products that can be used with Provox voice prostheses

- **Provox HME:** Heat- and Moisture Exchanger partially restores lost nasal functions (heating, humidifying, breathing resistance).
- **Provox FreeHands HME:** A Heat and Moisture Exchanger combined with an automatic tracheostoma valve, designed to allow hands-free speech.
- **Provox LaryTube:** A silicone tracheal cannula intended to support shrinking stomas while retaining other rehabilitation devices of the Provox System.
- **Provox LaryButton:** A soft, self-retaining silicone tracheal cannula intended to support shrinking stomas while retaining other rehabilitation devices of the Provox System.

4. Adverse events/Trouble shooting information

The following are possible risks associated with prosthetic voice rehabilitation:

Prosthesis related

Leakage through the valve – Leakage through the prosthesis may occur in association with Candida overgrowth of the valve flap or dislodgement of the radio opaque ring (valve seat). It may also occur if the esophageal flange of the prosthesis is not fully developed in the lumen of the esophagus or if the prosthesis is too short. Endoscopic examination should be performed to determine whether any of these phenomena has occurred. If the prosthesis is in the correct position and of the right length, and if cleaning of the prosthesis does not resolve leakage through the prosthesis, replacement of the prosthesis is indicated.

Candida overgrowth – Candida overgrowth of the prosthesis occurs in almost all patients. Radiotherapy, salivary content and dietary habits may influence how rapidly Candida can affect the silicone material. Treatment with antifungal agents may be considered as a preventive measure if Candida overgrowth causes leakage through the prosthesis or other incompetence of the valve.

Compatibility with antifungal medicine – An unknown variety of chemical substances may immediately influence the material properties of the device or may be absorbed by the material and subsequently released at a later time. Therefore, the introduction of antifungal medicines or drugs directly to or within close proximity of the voice prosthesis should be carefully assessed.

Laboratory tests show only minor negative influence on the function of the Provox Vega voice prosthesis and components when using the following antifungal medications:

- Nystatin
- Fluconazole
- Mikonazol

Puncture related

Aspiration of the prosthesis – Accidental aspiration of the voice prosthesis or other components of the voice rehabilitation system may occur. Immediate symptoms may include gagging, coughing, choking or wheezing. As with any other foreign body, complications from aspiration of a component may cause obstruction or infection. Complications may include pneumonia, atelectasis, bronchitis, lung abscess, bronchopulmonary fistula and asthma. If the patient can breathe, coughing may remove the foreign body. Partial airway obstruction or complete airway obstruction requires immediate intervention for removal of the object.

Ingestion of the prosthesis – Accidental ingestion of the voice prosthesis, or other components of the voice rehabilitation system, may occur. As with any other foreign body, the symptoms caused by ingestion of the prosthesis or component of the voice rehabilitation system depends largely on size, location, degree of obstruction (if any) and the length of time it has been present. Ingested components that have remained in the lower esophagus may be removed by esophagoscopy or observed for a short period of time. The object may pass spontaneously into the stomach; foreign bodies that pass into the stomach usually pass through the intestinal tract. Surgical removal of foreign bodies in the intestinal tract must be considered when bowel obstruction occurs, bleeding is present, perforation occurs or the object fails to pass through the intestinal tract.

Hemorrhage/Bleeding of the puncture – Slight bleeding from the edges of the TE-puncture may occur during replacement of the prosthesis and generally resolves spontaneously. Patients on anti-coagulant therapy, however, should be carefully evaluated for the risk of hemorrhage prior to placement or replacement of the prosthesis.

Infection and/or edema of the TE-puncture – Infection, granulation formation and/or edema of the puncture (e.g during radio therapy) may increase the length of the puncture tract. This may cause the prosthesis to be drawn inward and under the tracheal or esophageal mucosa. Or, inflammation or overgrowth of the esophageal mucosa may cause the prosthesis to protrude from the puncture. Temporary replacement of the prosthesis by a prosthesis with a longer shaft is then advisable. Treatment with broad-spectrum antibiotics with or without corticosteroids may be considered for treatment of the infection. If the infection does not resolve with antibiotics and/or corticosteroid intervention in the presence of the prosthesis, the prosthesis should be removed. In some cases stenting the puncture with a catheter might be considered. If the puncture closes spontaneously secondary to removal of the prosthesis, repuncture for insertion of a new prosthesis may be required.

Granulation around the puncture – Formation of granulation tissue around the TE-puncture has been reported at an incidence of 5%. Electrical, chemical, or laser cauterization of the area of granulation may be considered.

Hypertrophic scarring around the puncture – Bulging of the tracheal mucosa over the tracheal flange may occur if the prosthesis is relatively short. This excess tissue may be removed by using a laser (CO₂, or NdYAG). Alternatively, a prosthesis with a longer shaft can be used.

Protrusion/extrusion of the prosthesis – Protrusion of the prosthesis and subsequent spontaneous extrusion is sometimes observed during infection of the TE-puncture. Removal of the prosthesis is required to avoid dislodgement into the trachea. The puncture may close spontaneously secondary to the removal of the prosthesis. Repuncture may be necessary for insertion of a new prosthesis.

Leakage around the prosthesis – Transient leakage around the prosthesis may occur and may improve spontaneously. The most common reason is that the prosthesis is too long, which is solved by inserting a shorter prosthesis. A prosthesis of another diameter (upsizing) can also solve the leakage. It can also occur when the enlarged flange of the Provox Vega XtraSeal has not deployed fully into the esophageal lumen. No part of the enlarged esophageal flange should be visible alongside the shaft of the voice prosthesis and when rotated the voice prosthesis should move freely. Rotating the voice prosthesis while at the same time slightly pushing it towards the esophagus can help to unfold the enlarged flange completely. In case of uncertainty, use flexible endoscopy to confirm proper placement inside the esophagus.

If leakage does not improve (which it often does spontaneously), a Provox XtraFlange can be used. Other methods to solve leakage around the voice prosthesis are temporary removal of the prosthesis and the insertion of a cuffed tracheal cannula and/or nasogastric feeding tube to permit shrinkage of the puncture. Alternatively, a purse string suture with 3x0 absorbable material could be submucosally applied around the TE-puncture after removal of the prosthesis. The new prosthesis should be inserted and the suture should be tightened gently, enclosing the puncture walls around the prosthesis. If leakage around the prosthesis is intractable, more conservative measures, surgical closure of the puncture and subsequent re-puncture may be necessary.

Tissue damage – If the prosthesis is too short, too long, or is pushed frequently against the esophageal wall by a tracheal cannula, stoma button, or the patients' finger, damage of the puncture, tracheal and/or esophageal tissues may occur. In patients undergoing radiotherapy this may happen more easily. Inspect the conditions regularly to avoid severe damage.

5. Additional information

5.1 Compatibility with MRI, X-ray and radiation therapy

Provox voice prostheses have been tested and found to be compatible, except Provox ActiValve, with Magnetic Resonance Imaging (tested up to 3.0 Tesla), x-ray and radiation therapy (tested up to 70 Gy). The prosthesis can be left in the TE-puncture during the examination/radiation therapy session. Note that the radio-opaque valve seat is visible on x-ray.

5.2 Device lifetime

Depending on individual biological circumstances the device life varies, and it is not possible to predict the integrity of the device over a longer period of time. The material of the device will be affected by e.g. bacteria and yeast, and the structural integrity of the device will eventually deteriorate.

5.3 User assistance information

For additional help or information, please see back cover of the Instructions for Use for contact information.

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© Atos Medical AB, 2015.