



ADVANCED  
WOUND CARE

# OptiLock<sup>®</sup>

NON-ADHESIVE  
SUPER ABSORBENT  
WOUND DRESSING



Lock away exudate.

# The optimal dressing for heavily exuding wounds.

OptiLock wound dressing absorbs a remarkable amount of wound fluid—even under compression—and locks it away to help prevent maceration. A non-adherent facing minimizes disruption to the wound and allows for gentle removal that helps keep patients comfortable.

## Gentle wound care that's tough on exudate.



### Absorbs

- Super absorbent polymer core efficiently absorbs high volumes of exudate
- **Gently & efficiently absorbs,** even on the most heavily draining wounds



### Locks

- Inside the SAP core, exudate is converted into a gel that's kept firmly away from the wound, even under compression
- **Gently protects** against maceration, including the periwound area




### Adjusts

- Continually adapts the level of absorption to the amount of wound drainage
- **Gently provides** an optimal moist wound healing environment



### Protects

- Non-adherent wound contact layer prevents the dressing from sticking to the wound or periwound area
- **Gently releases** from the wound during dressing changes




**Wound Contact Layer**  
Non-adherent,  
for gentle removal.

**"Reverse Funnel" Design**  
Unique design of contact layer helps direct fluid in one direction—away from the wound.



# OptiLock dressing absorbs and retains more fluid.



## Dispersion Layer

Distributes and quickly transfers exudate evenly into the super absorbent core.

## Super Absorbent Core


Absorbs large amounts of exudate, which is converted to a gel and locked away from the wound.

## Blue Protective Backing

Creates a barrier, helping prevent strikethrough.


In manufacturer's laboratory testing, **OptiLock dressings** outperformed a leading silicone-faced foam dressing. OptiLock's remarkable absorption and fluid retention capabilities raise the bar for gentle wound care.

### Absorption Capacity Under Compression<sup>1</sup>



OptiLock **absorbs approximately three times more fluid** than a leading silicone-faced foam dressing.

### Fluid Loss Under Compression<sup>1</sup>



OptiLock dressings **retain substantially more fluid under compression** than a leading silicone-faced foam dressing.

# Our educational packaging takes the guesswork out of bandage application.

## OptiLock Educational Packaging

OptiLock features Medline's Educational Packaging (EP). The package focuses on patient safety and correct product application. Each package serves as a 2-Minute Course on Wound Care.™

## Uses

A primary dressing on partial-thickness and full-thickness wounds or a secondary dressing on full-thickness wounds, including:

- Ulcers (pressure, venous, neuropathic/diabetic)
- Postoperative wounds
- First- and second-degree burns

## Change Frequency

Can remain in place up to seven days, depending on the level of exudate.

## Ordering Information

Item No.	Description	Pkg.
MSC6433EP	3" x 3" (7.62 cm x 7.62 cm)	100/cs, 10/bx, 1 ea
MSC6444EP	4" x 4" (10.2 cm x 10.2 cm)	100/cs, 10/bx, 1 ea
MSC6455EP	5" x 5½" (12.7 cm x 14 cm)	100/cs, 10/bx, 1 ea
MSC64610EP	6½" x 10" (16.5 cm x 25.4 cm)	50/cs, 10/bx, 1 ea
MSC64812EP	8" x 12" (20.3 cm x 30.5 cm)	80/cs, 10/bx, 1 ea

To order by the box, add a "Z" to the end of the item number.  
To order by the each, add an "H" to the end of the item number.

### OPTILOCK

**NON-ADHESIVE SUPER ABSORBENT WOUND DRESSING**

OptiLock absorbs moderate to heavy drainage. Locks in fluid and protects skin around the wound. Gentle non-adherent wound contact layer.

**Usos:**  
Absorbe y retiene material drenado  
Para heridas con exposición parcial y total de tejido

**Frecuencia de cambio:**  
Dura hasta 7 días

1. Remove OptiLock from the package.

2. Place white side of dressing on the wound. Use a large enough dressing to cover the wound and surrounding skin.

3. Secure the dressing with tape or other secondary dressing such as elastic net.

**NOTE:** Do not cut OptiLock.

**CONTRAINDICATIONS:** Third degree burns. Known sensitivity to the product itself or its components. Do not use on non-exudating or minimally exudating wounds.

**STORAGE:** Protect from freezing. Avoid excessive heat. Sterile in unopened, undamaged package. Single use only.

1-800-MEDLINE (1-800-633-5463) www.medline.com/woundcare  
Manufactured in Austria for Medline Industries, Inc.,  
Mundelein, IL 60060 USA R010LRU

REORDER MSC644EP

**MEDLINE**

LOT

(01)30884389111271

10.2 cm x 10.2 cm DRESSING 10 PER BOX LATEX FREE STERILE  
8.4 cm x 8.9 cm PAD

**ABSORBENT POLYMER**

The product box and individual pouches feature step-by-step directions on how to use the dressing properly, making it easier to do the right thing.



**For exceptional wound care products** contact your Medline Representative or call 1-800-MEDLINE.



**Medline Industries, Inc.** Three Lakes Drive, Northfield, IL 60093  
Medline United States  
1-800-MEDLINE (633-5463)  
medline.com | info@medline.com

**Medline Canada**  
1-800-396-6996 | medline.ca | canada@medline.com  
**Medline México**  
01-800-831-0898 | medlinemexico.com | mexico@medline.com

Follow us