# maternity support soutien prénatal soporte de maternidad

# Loving Comfort

To Wear - façon de le porter - instrucciones di uso

## 1. Putting On The Abdominal Support Pad

Place the Abdominal Support Pad (ASP) in a lifting position under the abdomen. Fasten it in place by wrapping its strap around your body and touching the Velcro hook to the ASP.

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## 2. Position The ASP

Place the ASP as low under your abdomen as possible but not so low as to interfere with your legs when you walk or sit. Slide the ASP up slightly until it doesn't interfere with your legs when lifted as shown. The ASP is now positioned correctly.

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## 3. Attaching The Belt

Grip the belt securely at both ends as shown. Stretch it forward and touch it to the ASP. You can adjust the tension and support you feel according to the amount you stretch the belt. To prevent skin irritation, be sure the Velcro hook on the ends of the belt is entirely on the ASP and doesn't touch your skin.

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
## 4. Fasten The Upper Tummy Strap

For an extra sense of security, lightly stretch the tummy strap over and touch it to the Velcro loop on the opposite side. The tummy strap is not essential but many prefer to have it. You can also use the strap to help hold the back support up while you are positioning it. The choice to use it or remove it is yours.

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## 5. Final Step (Optional)

Lay on your back with knees drawn up. Release the ends of the belt, stretch it again, re-attach it to the ASP and stand up. This will lift the abdomen higher to achieve even greater relief from sagging weight and leg pain. Repeat this step each time you put on your support.

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## 6. Grows With Baby

Your Loving Comfort® Maternity Support is now in place and you're ready to go! As you and your baby grow, Loving Comfort® Maternity Support expands to accommodate by simply fastening further away from the center of the ASP.

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---

#### Washing instructions

hand wash in cold or warm water with mild detergent. Hang dry.  
Take care to cover the hook & loop so it doesn't tear the fabric  
as you handle it. Do not machine wash or dry.

#### Sizing

The loving comfort <sup>®</sup>Maternity Support is made of a  
latex-free, whisper soft elastic. For proper fit and maximum  
effectiveness, please use the following size chart. Refer to  
your pre-pregnancy dress size.

Size: (pre-pregnancydress Size)	
Size	fitsdress
SM	3-6
Md	7-16
IG	17-20

**Always consult your physician whenever you have questions about your health or your baby's health.**