12. Carefully examine cover before each use. Discard the pad if cover shows any sign of deterioration.
11. Do not use the cord as a handle.
10. Never pull this pad by the supply cord.
9. Do not sit on or crush pad – avoid sharp folds.
8. Do not use pins or other metallic means to fasten this pad in place.
7. The power unit is intended to be oriented in a vertical or floor mount position.

6. Burns can occur regardless of control setting: check skin under pad frequently.
 5. Do not use on areas of insensitive skin.
 4. This pad is not to be used on or by any individuals with certain disabilities, a sleeping or unconscious person, or a person with poor blood circulation, unless carefully attended.
 3. Do not use while sleeping.
 2. Do not use while sleeping.
- 1. ALWAYS UNPLUG THE UNIT FROM THE ELECTRICAL OUTLET IMMEDIATELY AFTER USING.**

READ ALL INSTRUCTIONS BEFORE USE

DANGER – To reduce the risk of burns, electric shock and fire, this product must be used in accordance with the following instructions:

IMPORTANT SAFETY INSTRUCTIONS DO NOT DESTROY

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- Keep away from children.
- Do not use on any unexplained pain, swollen muscles, or following a serious muscular injury before consulting your physician.
- **Consult your physician should you have any questions regarding therapeutic use of this heating pad.**
- This heating pad should not be used in the following cases unless directed by a physician: if you have diabetes, tuberculosis, benign or malignant tumors, phlebitis and thrombosis, hemorrhages, inflammations of the skin, undetermined calf pain or joint pain, or if you have poor circulation. Consult a physician before using if you have any kind of fungal infection, swelling, fracture or persistent pain.

Cautions

SAVE THESE INSTRUCTIONS

20. Use only supplied adapter with this pad.
19. Place pad on top of and not under the part of the body needing heat.
18. Use this pad only for its intended use as described in this manual.
17. Keep the pad away from heated surfaces.
16. Loop cord loosely when storing: tight wrapping may damage cord and internal parts.
15. Do not use in an oxygen-enriched environment.
14. Do not use this pad with liniment, salve or ointment preparations that contain heat-producing ingredients, as skin burns could result. Use recommended by your physician for the relief of pain.
13. Do not tamper with this pad in any way. There are no user-serviceable parts. If for any reason this pad does not function satisfactorily, return to authorized service center at address on back.

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The massage motors will automatically shut off after approximately 20 minutes. To restart, first set the massage switch back to OFF and then set the switch to LO or HI again. Please note the heating function operates independently of the timer.

Timer Feature

minutes at one time.

benefit, use should not exceed 20

HI levels for vibration. For maximum switch to choose between LO and Turn on the switch for heat. Set the slide into a standard 120V electrical outlet.

To operate the Therma Luxe Massaging Heating Pad, first plug the power jack into the hand control and plug the cord

IB-11069

LIMITED FIVE YEAR WARRANTY

Conair will repair or replace (at our option) your unit free of charge for 60 months from the date of purchase if the appliance is defective in workmanship or materials.

To obtain service under this warranty, return the defective product to the service center listed on back, together with your purchase receipt and \$3.00 for postage and handling. California residents need only provide proof of purchase and should call 1-800-366-0937 for shipping instructions. In the absence of a purchase receipt, the warranty period shall be 60 months from the date of manufacture.

ANY IMPLIED WARRANTIES, OBLIGATIONS, OR LIABILITIES, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, SHALL BE LIMITED IN DURATION TO THE 60-MONTH DURATION OF THIS WRITTEN, LIMITED WARRANTY. Some states do not allow limitations on how long an implied warranty lasts, so the above limitations may not apply to you.

IN NO EVENT SHALL CONAIR BE LIABLE FOR ANY SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES FOR BREACH OF THIS OR ANY OTHER WARRANTY, EXPRESS OR IMPLIED, WHATSOEVER.

Some states do not allow the exclusion or limitation of special, incidental, or consequential damages, so the above limitation may not apply to you.

This warranty gives you specific legal rights, and you may also have other rights, which vary from state to state.

Vistenos en www.conair.com
 @2011 Conair Corporation
 East Windsor, NJ 08520 / Glendale, AZ 85307
 HP08T
 11PA125999

Centro de servicio
 Conair Corporation
 Service Department
 7475 N. Glen Harbor Blvd.
 Glendale, AZ 85307

Para registrar su producto, visitenos en:
www.conair.com/registration

Please register this product at:
www.conair.com/registration

Service Center
 Conair Corporation
 Service Department
 7475 North Glen Harbor Blvd.
 Glendale, AZ 85307

Visit us on the web at: www.conair.com
 ©2011 Conair Corporation
 East Windsor, NJ 08520 / Glendale, AZ 85307

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