

P.T.Rail[™] Floor Installation Options A floor mounted column that the P.T. Rail[™] wall plate (both hinged & fixed models) can mount to. Floor mounted and embedded wall masts available.

Features and Benefits

Floor Mast.

A floor mounted column used when wall installation of a P.T. Rail™ is not possible or desirable. Height adjustable from 31-36" / 79-91cm from floor to top rail. Available in chrome and white. For stainless steel floor mast, the rail height is 36" / 91cm from floor to top of rail.

Wall Embedded Floor Mast.

Floor mounted column installs behind the wall during construction. Rail height is 33" / 84 cm from floor to top rail. Zinc coated alloy steel.

Customized Solutions.

For larger facility orders, P.T. Rail™ wall plates can be designed to accommodate architectural structures such as pre-existing wall plates and hole configurations. Call for details.





White & Chrome Floor Mast Specifications:

Rail Height Range	31-36" / 79-91 cm (from floor to top of rail).
Mast Height	39" / 99cm.
Clamp Plate	3.5" x 15" / 9 x 38 cm.
Floor Plate	4.25" x 6.25" / 11 x 16 cm.
Finishes	Powder coated (white) or satin chrome plated steel.

Stainless Steel Floor Mast Specifications:

Fixed Rail Height	36" / 91cm (from floor to top of rail).
Mast Height	39" / 99cm.
Floor Plate	4.25" x 6.25" / 11 x 16 cm.
Finishes	Satin stainless steel.

Wall Embedded Mast Specifications:

Rail Height	31-36" / 79-91 cm (from floor to top of rail).
Floor Plate	6" x 3.5" / 15 x 9 cm
Finishes / Material	Zinc plated steel.

HealthCraft Products Inc., 2790 Fenton Road, Ottawa, Ontario, Canada K1T 3T7 Toll Free: (888) 619-9992 Tel: (613) 822-1885 Fax: (613) 822-1886 Web: www.healthcraftproducts.com E-mail: info@healthcraftproducts.com

