

Purpose

Skil-Care's Econo-Gel Pads are designed to relieve pressure and increase comfort for wheelchair users. High-density foam layers on top and bottom provide stability and cushioning. The gel-filled center inserts add pressure-relieving flotation.

Important Notice

Econo-Gel Pads have safety ties that must be secured to the wheelchair to prevent the cushion from slipping off the wheelchair seat.

Easy Care

The vinyl cover may be wiped clean using a damp cloth and spray-type cleanser or mild disinfectant. The polyester cover can be laundered and dried at temperatures that do not exceed 180°F (82°C). Do not use bleach.



Econo-Gel Pads - Standard Series

Re-Order Number	Cover	Size
751150	Vinyl w/Polyester	16" x 16" x 2"
751151	Vinyl w/Polyester	18" x 16" x 2"
751153	Vinyl w/Polyester	20" x 16" x 2"
751154	Vinyl Only	16" x 16" x 2"
751155	Vinyl Only	18" x 16" x 2"
751157	Vinyl Only	20" x 16" x 2"
751162	Vinyl w/Sheepskin	16" x 16" x 2"
751163	Vinyl w/Sheepskin	18" x 16" x 2"
751164	Vinyl w/Sheepskin	20" x 16" x 2"
751166	LSII Cover	16" x 16" x 2"
751167	LSII Cover	18" x 16" x 2"
751168	LSII Cover	20" x 16" x 2"




Econo-Gel Pads - Bimini Series

Re-Order Number	Cover	Size
751180	Bimini Blue	16" x 16" x 2"
751181	Bimini Blue	18" x 16" x 2"

Econo-Gel Pads - Thin-Line Series

Re-Order Number	Cover	Size
751170	Vinyl w/Polyester	16" x 16" x 1.5"
751171	Vinyl w/Polyester	18" x 16" x 1.5"
751172	Vinyl Only	16" x 16" x 1.5"
751173	Vinyl Only	18" x 16" x 1.5"
751174	Vinyl w/Polyester	20" x 16" x 1.5"



Turn product insert over for more of Skil-Care's Pressure Control Cushions, Wheelchair Accessories and Bed Positioning products.


Skil-Care is the leader in Pressure Control Cushions, Wheelchair Accessories and Bed Positioning products. Choose Skil-Care, it's the comfortable way of thinking.

Cushions - Pressure Relief

Gel-Foam X-Cushion

Item#	Description	UM
751026	Vinyl Only - 18"x16"x2.5"	Ea.
751027	Vinyl w/Low-Shear cover - 18"x16"x2.5"	Ea.

- Combines stability and pressure relief in one cushion.
- Water based gel provides pressure relief and cooling for added protection.
- Heat sealed incontinence proof cover for added infection control.



X-Gel Overlay

Item#	Description	UM
914430	18"x16"x1"	Ea.
781065	Optional Low Shear II Cover	Ea.

- Four quadrant gel topper, enhances stability and pressure relief for any cushion.
- Heat sealed incontinence proof cover for added infection control.



Foam-Air Cushion

Item#	Description	UM
914251	16"x16"x2.5"	Ea.
914250	18"x16"x2.5"	Ea.
914252	20"x16"x2.5"	Ea.

- Air-quilt surface inflates or deflates to meet resident's pressure reduction and comfort requirements.
- High density foam base provides stability and prevents bottoming out.
- Air pump and cover included.



Wheelchair Accessories - Upper/Lower Body

Adjustable Lateral Support

Item#	Description	UM
706050	Small - Width Adjusts from 8"-12"	Ea.
706055	Medium - Width Adjusts from 12"-18"	Ea.
706059	Large/Bariatric - Width Adjusts from 16"-24"	Ea.

- Minimizes left-right leaning in wheelchairs and most other chairs.
- Foam-padded sides adjust in or out for custom fit.



Wheelchair Footrest Extender

Item#	Description	UM
703286	Fits 16", 18"x1" W Wheelchair	Ea.
703288	Fits 16", 18"x2" W Wheelchair	Ea.
703294	Fits 20", 22", 24"x1" W Wheelchair	Ea.
703296	Fits 20", 22", 24"x2" W Wheelchair	Ea.

- Provides foot support and leg cushioning.
- Raises footrest for improved posture and protects leg against contact with wheelchair.
- Prevents feet from slipping behind footrest.



Swing-Away Foot Support

Item#	Description	UM
703471	Right - Fits 16"-20" W Wheelchair	Ea.
703470	Left - Fits 16"-20" W Wheelchair	Ea.
703472	Fits 16"-20" W Wheelchair	Pair
703476	Right - Fits 22"-30" W Wheelchair	Ea.
703475	Left - Fits 22"-30" W Wheelchair	Ea.
703477	Fits 22"-30" W Wheelchair	Pair


- Cushioned back and side panels prevent feet from slipping off footrest.
- Swings away to allow easy exit from the wheelchair
- Available for standard-size and bariatric wheelchairs.



30° Degree Wedge

Item#	Description	UM
554025	7"x12"x24"	Pair

- Provides optimal 30-degree angle of incline for relieving pressure on hipbone and coccyx
- Wipes cleans
- Can be disinfected for infection control.



Bed Positioning

Roll Control Bolster

Item#	Description	UM
556010	Single	Ea.
556012	Double	Ea.

- Provides secure restraint free positioning for side-lying residents.
- Prevents resident from rolling off the bed.
- Covered in durable, wipe-clean vinyl.



Vinyl Bed Rail Pads

Item#	Description	UM
401015	Thru-View - 60"x15"x2"	Pair
401045	Thru-View - 72"x15"x2"	Pair
401010	Regular (none thru-view) - 60"x15"x1"	Pair
401140	Regular (none thru-view) - 72"x15"x2"	Pair

- Protects and cushions resident from bed rails.
- Velcro® closures fasten pads securely to bed rails.
- Covered in durable reinforced, wipe-clean vinyl.



For any questions please contact our customer service department.

29 Wells Ave, Yonkers, NY 10701 • 914.963.2040 • 800.431.2972 • fax: 914.963.2567

Visit us online at www.skil-care.com to view our full line of Skil-Care products