



Visco-Elastic Foam



Elastic Straps which wrap easily over any wheelchair foam cushion



Re-Order Number	Descriptions	Size	Units
914441	Visco Cushion Topper w/ Low Shear II Cover	18" x 16" x 1"	Each

Purpose

This thin-line covered Visco topper helps in upgrading any foam cushion. The effective pressure relieving Visco foam cushion and water proof cover are a must for enhancing any low-level foam cushion.

Important Safety Feature

The Visco Cushion Topper has elastic straps underneath that wrap around any foam wheelchair cushion. Which prevents the cushion from sliding off the foam cushion.

Easy Care

The Visco Cushion Topper wipes clean with a damp cloth and spray-type cleanser. It may be disinfected with a non-chlorine disinfectant. Do not use bleach or launder. The cloth cover can be cleaned with a wipe-on or spray disinfectant. Do not use bleach or launder.

This product is guaranteed for 6 months.