


Enteral Feeding Gravity Sets

McKesson enteral feeding gravity sets with ENFit connectors are to be used when no enteral feeding pump is available. The feeding sets are intended for administration of enteral nutrition, fluids or medication. They are available in large or small bore.

Features and Benefits

- All sets include a pre-attached feeding tube transition connector (for feeding tubes without the ENFit connection system)
- All sets also include an ENFit connector which eliminates I.V. misconnection, improving patient safety
- The connector is compliant with the ENFit ISO 80369-3 design standard which supports the goal to improve patient safety by reducing the risk of misconnections
- Graduated measurements and translucent material allows for the monitoring of fluid level
- Each gravity bag has a pouch to potentially hold ice for formulas or medications requiring cooler temperatures
- Roller clamp adjusts the flow rate of feeding
- Not made with DEHP, non-sterile
- Not made with natural rubber latex



ORDERING INFORMATION

ENTERAL FEEDING GRAVITY SETS

MFR #	DESCRIPTION	ATTRIBUTES	SPECS	INNER DIAMETER	OUTER DIAMETER	FRENCH SIZE BASED ON OUTER DIAMETER	QTY
16-702500	Standard Bore	With ENFit Connector and Transition Connector, Snap Cap, Non-Sterile	1200 mL	3.0 mm	4.1 mm	12 Fr	1 bag, 30 bags/cs
16-702505	Large Bore	With ENFit Connector and Transition Connector, Snap Cap, Non-Sterile	1200 mL	4.8 mm	7.0 mm	20-22 Fr	1 bag, 30 bags/cs

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

McKesson Medical-Surgical
9954 Mayland Drive
Richmond, VA 23233
mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health
800.347.2456

Home Medical Equipment
888.822.8111

Primary Care
866.625.2679

McKesson Brands Support
800.777.4908