# **Infant Heel Warmer**

McKesson Infant Heel Warmers are designed to help minimize discomfort and increase blood flow prior to taking a blood sample. Easy to activate, they provide a one-step, convenient alternative to the use of hot towels, quickly heating up to 105 degrees Fahrenheit and producing uniform heat consistently. Product requires no additional preparation. Designed for single use. Can be stored at room temperature.

## **Features and Benefits**

- Evenly distributed heat warms up the heel area prior to drawing blood sample
- Helps improve blood draw from infants
- Less squeezing reduces possibility of bruising and promotes faster healing
- Medical grade adhesive helps holds the heel warmer in place
- Not made with natural rubber latex



### ORDERING INFORMATION

DESCRIPTION	MFR #	QTY
Infant Heel Warmer, with Adhesive Tab, 2.75" x 4"	204	25ea/bx, 4bx/cs
Infant Heel Warmer, with Adhesive Tab, 3.75" x 5.5"	59-45HW	25ea/bx, 4bx/cs

## Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

#### **Never Compromise on Care SM**

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them — making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

© 2015 McKesson Medical-Surgical Inc.