## **Water Soluble Laundry Bags**

McKesson water soluble laundry bags are designed for one-time handling of heavily contaminated linens. The bag and contents can be placed in washer without additional handling. The bags dissolve completely in water temperature of  $140^{\circ} F$  ( $60^{\circ} C$ ) or above.

### **Features and Benefits**

- · Hot water soluble bags
- For one-time handling of heavily contaminated linens
- Bags dissolve completely in water temperature of 140°F (60°C) or above
- Bags will not damage laundry, equipment, or sewer system when used properly
- Bag tie included
- Not made with natural rubber latex



### ORDERING INFORMATION

MFR#	DESCRIPTION	ATTRIBUTES	SIZE	QTY
03-1347	Laundry Bag, 35 Gallon, Water Soluble	Heavy Duty, Clear	28 in x 39 in, Thickness: 1 mil	25/packs, 100 packs/cs
03-641A	Laundry Bag, 25 Gallon, Water Soluble	Clear	26 in x 33 in, Thickness: 0.8 mil	25/packs, 100 packs/
03-648A	Laundry Bag, 45 Gallon, Water Soluble	Clear	36 in x 39 in, Thickness: 0.8 mil	25/packs, 100 packs/

#### Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit within 12 months from the purchase date. McKesson Brands' satisfaction guarantee is limited to products purchased directly from McKesson. If a McKesson Brands product is not purchased directly from McKesson, any return request or request for refund must be directed to the selling party.

#### **Never Compromise on Care SM**

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

© 2020 McKesson Medical-Surgical Inc.

# **M**SKESSON

McKesson Medical-Surgical 9954 Mayland Drive Richmond, VA 23233 mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health 800.347.2456

Home Medical Equipment 888.822.8111

Primary Care 866.625.2679

McKesson Brands Support 800.777.4908