

# Introducing the QuickChange Wrap

Not a brief. Not a foley catheter. Not a condom catheter.

Instead, a simple, comfortable solution that prevents CAUTIs, IADs and staff injuries.

## Simple as...



# Field-tested with phenomenal results

In a survey, 91% of frontline caregivers who have tried the QuickChange Wrap say it is better than current choices, which are not really choices at all. Virtually all respondents agreed that the wrap is easy to deploy and change. Almost all said they would recommend the product to others.

"Less invasive and less strain on staff and patients."



#### For decades, nothing new has come along to care for men with urinary incontinence.

In hospitals, 15%-25% of patients are placed on catheters, joined to a urine output bag. Staff spend countless hours positioning, cleaning and worrying about Foley catheters, which can cause urinary tract infections, the most common healthcareassociated infection and a condition not reimbursed by payers.

In long-term care, staff change briefs several times per shift, which means lifting and rolling patients, resulting in injuries to backs, knees and feet. Patients and bed linens are often soiled, resulting in incontinence-associated dermatitis (IAD) and other skin problems. Administrators must keep multiple sizes of briefs and pads in inventory. In post-acute care, many patients are transitioned to condom catheters. Finding the right size condom catheter is not easy, and staff can waste several trying them out. Condom catheters often leak onto the patient and bed linens and can be painful when removed. IAD is a common outcome.

#### **Everything has changed.**

The QuickChange Wrap is the only solution designed exclusively for male UI. It is adjustable to fit every patient and can be placed in seconds without moving or even waking the patient. One wrap absorbs half of the average daily urine output, reducing incidence of CAUTIS, urethral scarring and IAD. Nurses are free to care for patients.

### What is it?

- It's a soft, cloth-like, multi-layer, absorbent wrap.
- Has a circular penile seal that prevents leaks and allows for all penis widths.
- Has a premium closure/hook that adheres tightly and safely.
- Has a top layer next to the skin that is soft and porous so the urine is immediately wicked away to the absorbent middle layer, where urine is captured as a gel, dry to the skin.

### Life just got better for:



PATIENTS: Fewer infections. Restored mobility. Less invasive.



STAFF: Easier than catheters. Fewer soiled bed linens. Places on patient in seconds.



ADMINISTRATORS: Reduce CAUTI occurences. Lower facility costs. Greater staff productivity.

#### Finally, a chance to replace these...with this:







For more information call 800.206.2816 or visit us at www.uimed.com