

EQUOS[®] wound care

For all your treatment needs

Five-layer silicone foam dressings

Secure yet gentle adherence to skin, with a foam core that ensures reliable absorption in wounds.

Item No.	Description	Pkg.
EQX2333	Five-layer silicone foam, 3" x 3" (7.5 x 7.5 cm)	10/bx, 100/cs
EQX2344	Five-layer silicone foam, 4" x 4" (10 x 10 cm)	10/bx, 100/cs
EQX2366	Five-layer silicone foam, 6" x 6" (15.2 x 15.2 cm)	10/bx, 100/cs
EQX2377	Five-layer silicone foam, Sacrum, 7" x 7" (17.8 x 17.8 cm)	5/bx, 40/cs
EQX2399	Five-layer silicone foam, Sacrum, 9" x 9" (22.8 x 22.8 cm)	5/bx, 25/cs



EQX2333

Silver Alginate wound dressings

Absorbent alginate primary dressing for partial-and full-thickness wounds.

Item No.	Description	Pkg.
EQX9012	Silver Alginate, 12" ribbon (30 cm)	5/bx, 50/cs
EQX9022	Silver Alginate, 2" x 2" (5 x 5 cm)	10/bx, 100/cs
EQX9045	Silver Alginate, 4" x 5" (10 x 12.7 cm)	10/bx, 100/cs
EQX9048	Silver Alginate, 4" x 8" (10 x 20.3 cm)	5/bx, 50/cs
EQX9066	Silver Alginate, 6" x 6" (15 x 15 cm)	5/bx, 50/cs



EQX9045

Alginate wound dressings

Absorbent alginate primary dressing for venous stasis ulcers, pressure ulcers, arterial ulcers, and diabetic ulcers.

Item No.	Description	Pkg.
EQX7012	Alginate, 12" ribbon (30 cm)	10/bx, 100/cs
EQX7022	Alginate, 2" x 2" (5 x 5 cm)	10/bx, 100/cs
EQX7044	Alginate, 4" x 4" (10 x 10 cm)	10/bx, 100/cs



EQX7022

Multi-layer compression systems

Delivers therapeutic compressions for the management of venous disease and associated edema.

Item No.	Description	Pkg.
EQX4300	Compression, 3-layer	1/bx, 12/cs
EQX4400	Compression, 4-layer	1/bx, 12/cs



EQX4300

Bordered gauze dressings

Absorbent island dressing with a soaker pad and non-woven fabric with an adhesive backing for use as a primary or secondary dressing.

Item No.	Description	Pkg.
EQX3223	Bordered gauze, 2" x 3.5" (5 x 9 cm)	50/bx, 500/cs
EQX3244	Bordered gauze, 4" x 4" (10 x 10 cm)	25/bx, 500/cs
EQX3266	Bordered gauze, 6" x 6" (15 x 15 cm)	25/bx, 500/cs



EQX3244

Wound cleanser

Formulated with a gentle surfactant to allow fast, thorough and gentle cleansing of dermal wounds. Spray delivery loosens wound debris.

Item No.	Description	Pkg.
EQX6008	EQUOS wound cleanser, 8oz. bottle with trigger sprayer	6/cs



EQX6008

Hydrogel dressing

Gel-based wound care dressing donates moisture to help create an optimal wound healing environment.

Item No.	Description	Pkg.
EQX6102	Hydrogel, 1-oz. tube (28g)	30/cs



EQX6102

Learn more. Contact your Medline Representative or visit medline.com.

